

contamination. The second clean-up step was not found necessary in the case of the Canadian blood samples used initially to develop the analytical procedure.

### 3.3.2.3 GC/MS Analysis

The cleaned-up extracts were chromatographed (without being derivatized) on a DB-5 capillary column, the eluted analyte being ionized in an oxygen atmosphere using a Townsend Discharge Ionizer. Minor fragmentation of the ions occurred, with both HT-2 and T-2 producing one high intensity ion ( $M + 32$ ) at 456 and 498  $m/z$  respectively, and lower intensity ions ( $M$ ) at 424 and 466  $m/z$  respectively. A Multiple Ion Detection (MID) profile of Monoacetoxyscirpenol (MAS), Diacetoxyscirpenol (DAS), T-2 toxin and HT-2 toxin was devised to monitor these four characteristic ions. Instrumental conditions for the gas chromatograph and mass spectrometer have already been given in section 3.3.1.

Each residue after clean-up was taken up in 50  $\mu$ l EtAc and 1  $\mu$ l was taken into a gas chromatographic syringe which had already been loaded with 1  $\mu$ l EtAc followed by a small plug of air. The sample was injected directly on-column using an on-column injector system (J & W Scientific Inc., Brockville, Ontario) and the column temperature was programmed. When data acquisition was completed, the GC column temperature was increased to remove contaminants, and 5  $\mu$ l hexane was injected to aid this process. GC column