


THE FORE-ARM LEVER WALKING STICK CRUTCH (PATENT)



CONSISTS of an almost horizontal handle, which is grasped in the palm of the hand, and of a slightly sloped spring which supports the forearm from the wrist to the elbow, where the arm hangs naturally. The greater part of the weight of the body bears on the hand and the wrist. Its use compels the active use of hands, arm and shoulder, and is thus an exercise for the upper limb. It brings back the joy of walking to the patient, who can easily and naturally use this crutch after a little practice, and half-a-dozen miles or more may be covered without fatigue. Stairs, 'buses and trams form no obstacles for the patient equipped with this invention.


DESCRIPTIVE LEAFLET SENT FREE.

Two Qualities,
18/6 and 28/6 per pair.

S. MAW, SON & SONS, LTD.,
(G.P.O. BOX No. 54),
7-12, ALDERSGATE STREET, LONDON, E.C. 1.

ALLEN & HANBURYS LIMITED, LONDON. AERO OPERATING URETHROSCOPE.

Wyndham Powell's Latest Model.



Made under the supervision of Mr. Wyndham Powell.

Descriptive Circular on application.

ALLEN & HANBURYS, LTD., Manufacturers of Surgical Instruments and Appliances,
48, WIGMORE STREET, LONDON, W. 1.
Works: Bethnal Green, E. 2.

CHAS. HEARSON & CO.'S SPECIALITIES.

Apparatus for Bacteriological, Pathological, Physiological, Chemical, Venereal, and Public Health Laboratories.

New Catalogue, 1918, printed in Five Languages, forwarded on application.

FACTORY AND SHOWROOMS—

68, WILLOW WALK, BERMONDSEY, LONDON, S.E. 1.

Bacteriological Incubator.

Cool Incubator.

CHARLES HEARSON & Co., Ltd., 235, Regent St., W. 1.

DUROGLASS.

STAR A BRAND

High Quality
Resistance Glassware
for all
BACTERIOLOGICAL
PURPOSES.

TWICE ANNEALED

BAIRD & TATLOCK (London), Ltd.,
14, Cross Street, Hatton Garden, E.C.



THE Surgical Tools

MANUFACTURED BY

MAYER & MELTZER

have a reputation for quality extending over
half-a-century.

Catalogues on request.

MAYER & MELTZER,

71, GT. PORTLAND STREET, LONDON, W. 1.

Telegrams: "TREPINE, WESDO."

Telephone: GERRARD 5574.