

For Diagnosis of Typhold Fever


Our Typhoid Agglutometer simplifies the Widal test, obviating the use of themicroscope and the fresh live culture of typhoid bacilli. It is

fully equal to the former method in delicity and accuracy, and much ruperior in convenience.

The apparatus consists
of four tubes of a sterile
permanent suspension of
dead typhoid bacilli, a tubo

of blood-serum diluent, an empty tube for collecting the blood, a needle for puncture, and a pipette for withdrawing the separated serum. It contains material for fifteen to thirty tests.

Directions Accompany Each Package Circular Free



PARKE DAVIS & COMPANY

LABORATORIES: DETROIT, MICH., U.E.A.; WALKERVILLE, ONT.; HOUNBLOW, ENG.

Branches: New York, Chicago, at. Louis, Boaton, Baltimore, New Orleans, Kansas City, Indianapolis, Minneapolis, Memphis; London, End'; Montheall, Que; Bydney, N. S. W.; 37, Pettresurg, Russis; Simla, Hodal Toxid, Japan.