



LONGITUDINAL SECTION OF THE CYSTOSCOPE.

FIG. 1.— With the obturator in place the sharp angle near the tip of the instrument is completely removed. The lamp is large and is easily withdrawn from the tip of the instrument. There are no wires to get out of order. The electrical attachment is indicated and may be made by an interlocking device if so desired. The wires *a* and *b*, if so desired, can be covered by rubber tubing, allowing of their sterilization. At Fig. 2 in *a* and *b* the other ends of the electric wires which are usually six to eight feet long are shown.