

GLYCERIN

Inspector's Report (Is not an expression of opinion).	Results of Analysis						Number of Sample	Remarks and Opinion of the Chief Analyst
	Refractive Index of Sample	Refractive Index of Distilled water	Difference	% Glycerol (Lens Tables)	Specific Gravity of sample at 20° C.	% Glycerol from Sp. Gr.		
F. X. W. BELAND, INSPECTOR								
.....	1.4700	1.3340	0.1360	95.5			56859	Genuine
.....	1.4710	1.3340	0.1370	96.3			56860	Genuine
.....	1.4727	1.3340	0.1387	97.5			56861	Genuine
.....	1.4730	1.3335	0.1395	97.8			56862	Genuine
.....	1.4635	1.3330	0.1307	92.0	1.2435	93.0	56863	Doubtful
.....	1.4550	1.3330	0.1220	86.0	1.2296	87.6	56864	Adulterated by dilution with water.
.....	1.4688	1.3335	0.1353	95.0			56865	Genuine
.....	1.4730	1.3340	0.1390	97.5			56866	Genuine
.....	1.4745	1.3340	0.1405	98.6			56867	Genuine
.....	1.4741	1.3335	0.1406	98.6			56868	Genuine

A. PELLETIER, T. GENDREAU, TEMPORARY INSPECTORS.

.....	1.4660	1.3330	0.1330	93.5			4366	Doubtful.
.....	1.4672	1.3330	0.1342	94.3			4368	Genuine
.....	1.4700	1.3330	0.1370	96.3			4372	Genuine
.....	1.4700	1.3330	0.1370	96.3			4380	Genuine
.....	1.4715	1.3330	0.1385	97.3			4386	Genuine
.....	1.4692	1.3330	0.1362	95.6			4387	Genuine
.....	1.4688	1.3330	0.1358	95.4			4399	Genuine
.....	1.4440	1.3330	0.1110	78.4	1.2347	Invert sugar 10.08 Sucrose 2.91 Total 12.99	5764	Adulterated by addition of sugar syrup
.....	1.4700	1.3330	0.1370	96.3			5765	Genuine
.....	1.4702	1.3330	0.1372	96.3			5766	Genuine
.....	1.4479	1.3330	0.1149	81.0	1.2474	Invert Sugar 10.92 Sucrose 5.59 Total 16.51	5775	Adulterated by addition of sugar syrup
.....	1.4711	1.3330	0.1381	97.0			5778	Genuine
.....	1.4705	1.3330	0.1375	96.7			5779	Genuine
.....	1.4460	1.3330	0.1130	79.8	1.2137	81.7	5780	Adulterated by added water
.....	1.4550	1.3330	0.1220	86.0	1.2288	87.3	4902	Adulterated by added water
.....	1.4445	1.3330	0.1115	78.8	1.2128	81.40	4996	Adulterated by added water
.....	1.4690	1.3330	0.1360	94.9			4999	Genuine
.....	1.4721	1.3330	0.1391	97.6			5002	Genuine
.....	1.4700	1.3330	0.1371	96.3			5006	Genuine
.....	1.4721	1.3330	0.1391	97.6			5009	Genuine
.....	1.4690	1.3330	0.1360	95.5			5011	Genuine
.....	1.4682	1.3330	0.1352	95.0			7302	Genuine
.....	1.4725	1.3330	0.1395	98.0			7307	Genuine