


What Clinical Trial Shows

- 1. Duration of the disease is cut short.
- 2. The most marked and disagreeable symptoms are ameliorated.
- 3. Temperature is reduced.
- 4. Delirium is terminated or much abated.
- 5. Tympanites is almost wholly avoided.
- 6. Diarrhea is checked.

remarkable antisepticpowerful, in that it is capable of dostroying any disease-germ with

which it can be brought into contact; safe, in that it may be taken inter-

nally, in aqueous solution, in quantities limited



Full Literature Sent Free on Request



For Diagnosis of Typhoid Fever

Our Typhoid Agglutometer simplifies the Widal test, obviating the use of the microscope and the fresh live culture of typhoid bacilli. It is

fully equal to the former method in delicacy and accuracy, and much superior in convenience.

> The apparatus consists of four tubes of a sterile permanent suspension of dead typhoid bacilli, a tube

of blood-serum diluent, an empty tube for collecting the blood, a needle for puncture,

and a pipette for withdrawing the separated serum. It contains material for fifteen to thirty tests.

Directions Accompany Each Package Circular Free



W YORK, CHICAGO, ST. LOUIS, BOSTON, BALTIMORE, NEW OFILEANS, KANSAS CITY, INDIANAPOLIS, MINNEAPOLIS, MEMPHIS; LONDON, ENG.; MONTREAL, QUE.; SYDNEY, N.S.W.; ST. PETERSBURG, RUSSIA; SIMLA, INDIA; TOKIO, JAPAN.