

Actien-Gesellschaft für Anilin-Fabrikation, Berlin S.O.

Photographic-Department.

TRADE



MARK

## The Developers:

**Rodinal :**

Highly concentrated developing solution for all kinds of photographic work. Diluted with 25 to 30 times its volume of water, it is ready for use.

**Eikonogen :**

Working very harmoniously, specially adapted for portrait and instantaneous photography.

**Amidol :**

Needing no alkali, can be well applied to all kinds of photographic work.

**Diogen :**

Elastic ! Powerful ! Handy ! Therefore suited for cases, in which it is doubtful, whether the plate has been correctly exposed.

**Diphenal :**

Concentrated solution particularly adapted for developing freely exposed negatives of landscapes.

**GLYCIN.**

**METOL.**

**ORTOL.**

**EIKONOGEN**—Cartridges or glass-tubes.

**AMIDOL**—Cartridges or glass-tubes.

Sold by all dealers. Ask your dealer for formulae.