




THREE POINTS OF Superiority

Combining other features of excellence IN THE Stearns "SIMPLEX" Syringe

- 
1. No glass to break. Syringe is sealed at the stem by a sterile cap. Slip off this cap, attach the needle and

you are ready to inject the serum. No other antitoxin syringe than ours has this feature.

2. A heavy rubber washer protects and strengthens the flange, and gives the fingers a firmer grip, making this an absolutely safe "one-hand" syringe. No other syringe has this feature.
3. Piston is firmly fixed to piston-rod, and under complete control of operator. No other syringe has this feature.

The "Simplex" is the climax of, simplicity in use as well as in construction. It is the only completely successful piston syringe on the market.

Our "Simplex" is a modification of the Prescott syringe, patented and assigned to us in 1901. Prescott's was the first piston syringe serum container ever made or patented.

Stearns' Diphtheritic Antitoxin is manufactured under Federal authority and inspection, and we believe it to be unequalled for purity and activity. Our laboratories and stables are models of cleanliness and scientific equipment, and our products are subjected to the most rigid bacteriologic and physiologic tests before being marketed.

When prescribing antitoxin tell the druggist "It must be Stearns'" and specify "Simplex" or "Syr-Bulb," according to your choice of syringes

FREDERICK STEARNS & CO.

DETROIT, MICH., U. S. A.

WINDSOR, ONT.

LONDON, ENG.

NEW YORK CITY



STEARNS QUALITY

FOUNDED
1855