

For Diagnosis of Typhoid Fever

Our Typhoid Agglutometer simplifies the Widel test, obviating the use of the microscope and the fresh live culture of typhoid bacilli. It is

> fully equal to the former method in delicacy and accuracy, and much superior in convenience.

> > The apparatus consists of four tubes of a sterile permanent suspension of

dead typhold bacilli, a tube

of blood-scrum diluent, an empty tube for collecting the blood, a needle for puncture,

and a pipette for withdrawing the separated serum. It contains material for fifteen to thirty tests.

Directions Accompany Each Package



MEMPHIS; LONDON, RNQ.; MONTREAL, QUE.; SYDNEY, N.S.W.; PETERBBURG, AUSGIA: SIMLA, INDIC! TOKIO, JAPAN.